


[Kd Validated] Anti-Echs1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30084

Item Number: CM0010178



Introduction

[KD Validated] ECHS1 Rabbit mAb for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat ECHS1. Recombinant immunogen; targets enoyl-CoA hydratase in fatty acid oxidation and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ECHS1 Rabbit mAb
Remarks/Alias	SCEH; mECH; mECH1; ECHS1D
Species	Human
GeneID (Human)	1892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-290 of human ECHS1 (NP_004083.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30084
Protein Weight	31kDa
Shipping	Ice bag