


[Kd Validated] Anti-Egfr Rabbit Polymab® For Wb, If/Icc, Elisa - P00533

Item Number: CM0005827



Introduction

Recombinant rabbit monoclonal EGFR antibody [KD Validated] for Western blot, IF/ICC, and ELISA. Reacts with human EGFR, a key receptor tyrosine kinase involved in cancer and signal transduction. Ideal for oncology research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] EGFR Rabbit PolymAb®
Remarks/Alias	EGFR; ERBB; ERBB1; HER1; NISBD2; PIG61; mENA; epidermal growth factor receptor
Species	Human
Gene ID (Human)	1956
Immunogen	Recombinant protein of human EGFR
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	44kDa/69kDa/77kDa/134kDa
Shipping	Ice bag