

[Kd Validated] Anti-Eif4E Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P06730

Item Number: CM0023617



Introduction

Rabbit polyclonal antibody against human eIF4E, a translation initiation factor. Validated for WB, IHC-P, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat eIF4E. Ideal for studying cap-dependent translation and mRNA export.

[Learn More](#)

Parameter	Specification
Product Name	[KD Validated] eIF4E Rabbit pAb
Remarks/Alias	CBP; EIF4F; AUTS19; EIF4E1; eIF-4E; EIF4EL1; 4E
Species	Human
GeneID (Human)	1977
GeneID	1977
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human eIF4E (NP_001959.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06730
Protein Weight	25kDa
Shipping	Ice bag