

[Kd Validated] Anti-Erk1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27361

Item Number: CM0013511



Introduction

[KD Validated] Rabbit polyclonal antibody against human ERK1, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for MAPK/ERK signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ERK1 Rabbit pAb
Remarks/Alias	ERK1, ERT2, ERK-1, PRKM3, P44ERK1, P44MAPK, HS44KDAP, HUMKER1A, p44-ERK1, p44-MAPK, K1
Species	Human
Gene ID	5595
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human ERK1 (NP_002737.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27361
Protein Weight	43kDa
Shipping	Ice bag