


[Kd Validated] Anti-Ferritin Light Chain (Ftl) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P02792

Item Number: CM0002011



Introduction

Rabbit monoclonal antibody targeting human ferritin light chain (FTL). Validated for WB, IHC-P, ELISA; cross-reacts with mouse, rat. Useful for iron homeostasis and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Ferritin Light Chain Rabbit mAb
Remarks/Alias	LFTD; NBIA3; Ferritin Light Chain
Species	Human
Gene ID (Human)	2512
Gene ID	2512
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human Ferritin Light Chain (NP_000137.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02792
Protein Weight	20kDa
Shipping	Ice bag