

[Kd Validated] Anti-Galc Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P54803

Item Number: CM0006996

Introduction

Rabbit monoclonal antibody against human galactocerebrosidase (GALC), validated for Western blot, immunohistochemistry, and ELISA. Detects human, mouse, and rat GALC (~77 kDa). Ideal for lysosomal storage disorder research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] GALC Rabbit mAb
Remarks/Alias	GALC; galactocerebrosidase
Species	Human
GeneID (human)	2581
GeneID	2581
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human GALC (NP_000144.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54803
Protein Weight	77kDa
Shipping	Ice bag