

[Kd Validated] Anti-Glucose-6-Phosphate Isomerase Rabbit Mab For Wb, Elisa - P06744

Item Number: CM0010683



Introduction

KD validated rabbit monoclonal antibody targeting human glucose-6-phosphate isomerase (GPI). Involved in glycolysis, and functions as cytokine AMF and neuroleukin NLK. Validated for Western blot (WB) and ELISA.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] glucose-6-phosphate isomerase Rabbit mAb
Remarks/Alias	AMF; NLK; PGI; PHI; GNPI; SA36; SA-36; [KD Validated] glucose-6-phosphate isomerase
Species	Human
Gene ID	2821 (Human)
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Glucose 6 phosphate isomerase (P06744)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06744
Protein Weight	63kDa
Shipping	Ice bag