

Kd Validated Anti-Gm130 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q08379

Item Number: CM0010632



Introduction

KD validated rabbit monoclonal antibody targeting GM130 (GOLGA2) for Western blot, immunofluorescence, IHC, and ELISA. Detects human, mouse, and rat. Suitable for Golgi apparatus and cell division research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] GM130 Rabbit mAb
Remarks/Alias	GM130; DEDHMB
Species	Human
Gene ID (Human)	2801
Gene ID	2801
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human GM130 (NP_004477.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08379
Protein Weight	113kDa
Shipping	Ice bag