

[Kd Validated] Anti-Gpld1 Rabbit Polyclonal Antibody For Wb, Elisa - P80108

Item Number: CM0020054

Introduction

High-quality [KD Validated] GPLD1 rabbit polyclonal antibody validated for WB and ELISA, recognizing human and mouse GPLD1 (~92kDa). Ideal for studying GPI anchor cleavage and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] GPLD1 Rabbit pAb
Remarks/Alias	PLD; GPIPLD; PIGPLD; GPIPLDM; PIGPLD1; [KD Validated] GPLD1
Species	Human
Gene ID (Human)	2822
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-160 of human GPLD1 (NP_001494.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P80108
Protein Weight	92kDa
Shipping	Ice bag