


[Kd Validated] Anti-Gpx4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36969

Item Number: CM0022677



Introduction

[KD Validated] rabbit polyclonal antibody against GPX4, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Essential for studying ferroptosis and oxidative stress. 22kDa protein.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] GPX4 Rabbit pAb
Remarks/Alias	MCSP; SMDS; GPx-4; PHGPx; snGPx; GSHPx-4; snPHGPx; GPX4
Species	Human
GeneID (Human)	2879
GeneID	2879
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 74-197 of GPX4 (NP_002076.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36969
Protein Weight	22kDa
Shipping	Ice bag