


[Kd Validated] Anti-Hsp90B1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P14625

Item Number: CM0008758



Introduction

KD-validated rabbit monoclonal antibody against HSP90B1 (GRP94), for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying ER chaperone function and Wnt/TLR signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] HSP90B1 Rabbit mAb
Remarks/Alias	ECGP; GP96; TRA1; GRP94; HEL35; HEL-S-125m
Species	Human
GeneID (Human)	7184
GeneID	7184
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 616-715 of human HSP90B1 (NP_003290.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14625
Protein Weight	92kDa
Shipping	Ice bag