

Kd Validated Anti-IL13 Monoclonal Antibody For Wb And Elisa - P35225

Item Number: CM0006451



Introduction

KD-validated rabbit monoclonal antibody against human IL13, suitable for Western blot and ELISA. Recombinant immunogen (aa 25-146). Recognizes human IL13. For allergic inflammation and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] IL13 Rabbit mAb
Remarks/Alias	P600; IL-13; [KD Validated] IL13
Species	Human
Gene ID	3596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-146 of human IL13. (NP_002179.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35225
Protein Weight	16kDa
Shipping	Ice bag