

[Kd Validated] Anti-Integrin Alpha 5 (Itga5/Cd49E) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P08648

Item Number: CM0007158



Introduction

KD-validated rabbit monoclonal antibody targeting human Integrin alpha-5 (ITGA5/CD49e). Confirmed for WB, IF/ICC, ELISA. Ideal for fibronectin receptor and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Integrin alpha 5 (ITGA5/CD49e) Rabbit mAb
Remarks/Alias	FNRA; CD49e; VLA-5; VLA5A; [KD Validated] Integrin alpha 5 (ITGA5/CD49e)
Species	Human
Gene ID	3678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-995 of human Integrin alpha 5 (ITGA5/CD49e) (NP_002196.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag