


[Kd Validated] Anti-Isg15 Polyclonal Antibody For Wb, Elisa - P05161

Item Number: CM0019410



Introduction

[KD Validated] Rabbit polyclonal antibody against human ISG15 (G1P2, UCRP), validated for WB and ELISA. Suitable for innate immunity, antiviral defense, and ISGylation research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ISG15 Rabbit pAb
Remarks/Alias	G1P2; IP17; UCRP; IFI15; IMD38; hUCRP; [KD Validated] ISG15
Species	Human
Gene ID (Human)	9636
Gene ID	9636
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human ISG15 (NP_005092.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P05161
Protein Weight	18kDa
Shipping	Ice bag