


[Kd Validated] Anti-Lc3B Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Gzq8

Item Number: CM0026319



Introduction

KD validated anti-LC3B rabbit polyclonal antibody from CamelBio. Designed for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with human, mouse, and rat LC3B. A core tool for studying autophagy and mitophagy pathways.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LC3B Rabbit pAb
Remarks/Alias	LC3B; ATG8F; MAP1LC3B-a; MAP1A/1BLC3; 3B
Species	Human
GeneID (Human)	81631
GeneID	81631
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LC3B (NP_073729.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZQ8
Protein Weight	15kDa
Shipping	Ice bag