

[Kd Validated] Anti-Lxr Beta/Ner Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P55055

Item Number: CM0007825



Introduction

Rabbit monoclonal antibody against human LXR beta/NER (NR1H2), validated for western blot, immunoprecipitation, and ELISA. Ideal for cholesterol metabolism and nuclear receptor research. Cross-reacts with human. - P55055

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] LXR beta/NER Rabbit mAb
Remarks/Alias	NER; UNR; LXRβ; LXR-b; NER-I; RIP15
Species	Human
GeneID (Human)	7376
GeneID	7376
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LXR beta/NER (NP_009052.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P55055
Protein Weight	51kDa
Shipping	Ice bag