

## [Kd Validated] Anti-Lypla1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O75608

Item Number: CM0010553



### Introduction

Rabbit monoclonal antibody against human LYPLA1 (APT1/LPL1), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets acyl-protein thioesterase 1 involved in protein depalmitoylation and lysophospholipid metabolism.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] LYPLA1 Rabbit mAb                                                                                  |
| Remarks/Alias    | APT1; LPL1; APT-1; LPL-I; hAPT1                                                                                   |
| Species          | Human                                                                                                             |
| Gene ID          | 10434                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human LYPLA1 (NP_006321.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | O75608                                                                                                            |
| Protein Weight   | 25kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |