


[Kd Validated] Anti-Mad1/Mad1L1 Polyclonal Antibody For Wb, Elisa - Q9Y6D9

Item Number: CM0023179



Introduction

Rabbit polyclonal antibody against human MAD1/MAD1L1, validated for Western blot and ELISA. Suitable for spindle assembly checkpoint and mitosis studies, with cross-reactivity to human samples.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MAD1/MAD1L1 Rabbit pAb
Remarks/Alias	MAD1; MVA7; PIG9; TP53I9; TXBP181; [KD Validated] MAD1/MAD1L1
Species	Human
GeneID (Human)	8379
GeneID	8379
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 351-718 of human MAD1/MAD1L1 (NP_003541.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y6D9
Protein Weight	83kDa
Shipping	Ice bag