


Kd Validated Anti-Mccc1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q96Rq3

Item Number: CM0004749



Introduction

KD validated rabbit monoclonal antibody against MCCC1, recognizing human, mouse, rat species. Suitable for WB, IF/ICC, IP, ELISA. Ideal for leucine catabolism and 3-methylcrotonylglycinuria research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MCCC1 Rabbit mAb
Remarks/Alias	MCCA; MCC-B; MCCCalpha; [KD Validated] MCCC1
Species	Human
GeneID (human)	56922
GeneID	56922
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RQ3
Protein Weight	80kDa
Shipping	Ice bag