


Kd Validated Anti-Me1 Rabbit Polyclonal Antibody For Wb And Elisa - P48163

Item Number: CM0023243



Introduction

KD validated rabbit polyclonal antibody targeting human ME1 (NADP-dependent malic enzyme) for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for metabolism and cancer research. SWISS: P48163.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ME1 Rabbit pAb
Remarks/Alias	MES; HUMNDME; [KD Validated] ME1
Species	Human
Gene ID (Human)	4199
Gene ID	4199
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 472-572 of human ME1 (NP_002386.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48163
Protein Weight	64kDa
Shipping	Ice bag