

[Kd Validated] Anti-Mettl3 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q86U44

Item Number: CM0002609



Introduction

Rabbit monoclonal antibody against human METTL3 (N6-adenosine-methyltransferase catalytic subunit), validated for WB, IP, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for epitranscriptomics research involving m6A RNA modification and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] METTL3 Rabbit mAb
Remarks/Alias	M6A; IME4; Spo8; MT-A70; hMETTL3; METTL3
Species	Human
Gene ID (Human)	56339
Gene ID	56339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-250 of human METTL3 (Q86U44).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86U44
Protein Weight	64kDa
Shipping	Ice bag