


[Kd Validated] Anti-Mif Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P14174

Item Number: CM0022834



Introduction

Rabbit polyclonal antibody against human MIF, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. For macrophage migration inhibitory factor research in inflammation and immunity.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MIF Rabbit pAb
Remarks/Alias	GIF; GLIF; MMIF; IF
Species	Human
GeneID (Human)	4282
GeneID	4282
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 2-115 of human MIF (NP_002406.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14174
Protein Weight	12kDa
Shipping	Ice bag