

[Kd Validated] Anti-Mrp4/Abcc4 Rabbit Monoclonal Antibody For Wb, Elisa - O15439

Item Number: CM0006994



Introduction

KD validated rabbit monoclonal antibody against human MRP4/ABCC4, validated for WB and ELISA. Exhibits cross-reactivity with mouse and rat. Ideal for multidrug resistance, transport studies, and cyclic nucleotide signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRP4/ABCC4 Rabbit mAb
Remarks/Alias	MRP4; MOATB; MOAT-B; C4
Species	Human
Gene ID (Human)	10257
Gene ID	10257
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1067-1325 of human MRP4/ABCC4 (NP_005836.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15439
Protein Weight	150kDa
Shipping	Ice bag