


[Kd Validated] Anti-Mrpl13 Rabbit Monoclonal Antibody (Wb, Ihc-P, Elisa) - Q9Byd1

Item Number: CM0005285



Introduction

KD validated MRPL13 rabbit monoclonal antibody recognizing human, mouse, rat MRPL13 protein. Suitable for WB, IHC-P, ELISA. Targets mitochondrial ribosomal protein uL13m (Q9BYD1). Ideal for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRPL13 Rabbit mAb
Remarks/Alias	L13; L13A; L13mt; RPL13; RPML13
Species	Human
GeneID (Human)	28998
GeneID	28998
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 70-178 of human MRPL13 (NP_054797.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYD1
Protein Weight	21kDa
Shipping	Ice bag