

Kd Validated Anti-Mtmr3 Polyclonal Antibody For Wb And Elisa - Q13615

Item Number: CM0015989

Introduction

KD Validated rabbit polyclonal antibody against human MTMR3 (ZFYVE10, FYVE-DSP1). Suitable for Western blot (WB) and ELISA. Targets a lipid phosphatase critical for PI3P metabolism and cytokinesis.

[Learn More](#)

Parameter	Details
Product Name	[KD Validated] MTMR3 Rabbit pAb
Remarks/Alias	ZFYVE10; FYVE-DSP1; [KD Validated] MTMR3
Species	Human
GeneID (Human)	8897
GeneID	8897
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 640-810 of human MTMR3 (NP_066576.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13615
Protein Weight	134kDa
Shipping	Ice bag