

[Kd Validated] Anti-Mtor Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P42345

Item Number: CM0006469



Introduction

KD validated mTOR rabbit monoclonal antibody for WB, IHC-P, ELISA. Reacts with human, mouse, rat mTOR; targets central metabolism and growth regulator. Ideal for research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] mTOR Rabbit mAb
Remarks/Alias	FRAP, FRAP1, FRAP2, RAFT1, RAPT1
Species	Human
Gene ID (Human)	2475
Gene ID	2475
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42345
Protein Weight	289 kDa
Shipping	Ice bag