

Kd Validated Anti-Ndufa13/Grim19 Polyclonal Antibody For Wb, Elisa - Q9P0J0

Item Number: CM0025457



Introduction

KD Validated rabbit polyclonal antibody against NDUFA13/GRIM19, validated for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for mitochondrial complex I and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] NDUFA13/GRIM19 Rabbit pAb
Remarks/Alias	B16.6; CDA016; CGI-39; GRIM19; GRIM-19; MC1DN28; [KD Validated] NDUFA13/GRIM19
Species	Human
GeneID (human)	51079
GeneID	51079
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human NDUFA13/GRIM19 (NP_057049.5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P0J0
Protein Weight	17kDa
Shipping	Ice bag