

## Kd Validated Anti-Nf2 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P35240

Item Number: CM0014198



### Introduction

Rabbit polyclonal antibody targeting human NF2 (Merlin), a tumor suppressor in the Hippo pathway. Validated for WB, IF/ICC, IP, and ELISA; cross-reacts with mouse and rat. Suitable for cancer and neurofibromatosis type 2 research.

[Learn More](#)

| Parameter        | Value                         |
|------------------|-------------------------------|
| Product Name     | [KD Validated] NF2 Rabbit pAb |
| Remarks/Alias    | ACN; SCH; BANF; merlin-1; F2  |
| Species          | Human                         |
| GeneID (Human)   | 4771                          |
| GeneID           | 4771                          |
| Immunogen        | Recombinant protein           |
| Source           | Rabbit                        |
| Category         | Polyclonal Antibodies         |
| Application      | WB, IF/ICC, IP, ELISA         |
| Cross Reactivity | Human, Mouse, Rat             |
| SWISS            | P35240                        |
| Protein Weight   | 70kDa                         |
| Shipping         | Ice bag                       |