

[Kd Validated] Anti-Nf2/Merlin Rabbit Monoclonal Antibody For Wb, Ihc-P And Elisa - P35240

Item Number: CM0004018



Introduction

KD validated recombinant rabbit monoclonal antibody against human NF2/Merlin, suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. For tumor suppressor and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] NF2/Merlin Rabbit mAb
Remarks/Alias	ACN; SCH; BANF; merlin-1; in
Species	Human
GeneID (Human)	4771
GeneID	4771
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 347-595 of human NF2/Merlin (NP_000259.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35240
Protein Weight	70kDa
Shipping	Ice bag