

## Kd Validated Anti-Pdxk Polyclonal Antibody For Wb, If/Icc, Elisa - O00764

Item Number: CM0021645



### Introduction

KD validated rabbit polyclonal antibody against human pyridoxal kinase (PDXK), suitable for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for vitamin B6 metabolism research.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated]PDXK Rabbit pAb                                                                        |
| Remarks/Alias    | PKH; PNK; HMSN6C; PRED79; C21orf97; HEL-S-1a; C21orf124; PDXK                                        |
| Species          | Human                                                                                                |
| GeneID (Human)   | 8566                                                                                                 |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 213-312 of human PDXK (NP_003672.1) |
| Source           | Rabbit                                                                                               |
| Category         | Polyclonal Antibodies                                                                                |
| Application      | WB, IF/ICC, ELISA                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                    |
| SWISS            | O00764                                                                                               |
| Protein Weight   | 35kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |