

[Kd Validated] Anti-Phf11 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Uil8

Item Number: CM0006400

Introduction

KD validated rabbit monoclonal antibody against PHF11 (BCAP), suitable for WB and ELISA. Detects PHF11 in human, mouse, and rat samples. Useful for studying Th1-type immune responses and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] PHF11 Rabbit mAb
Remarks/Alias	APY; BCAP; IGEL; IGER; IGHER; NYREN34; NY-REN-34; [KD Validated] PHF11
Species	Human
GeneID (Human)	51131
GeneID	51131
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 40-321 of human PHF11 NP_001035533.1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIL8
Protein Weight	38 kDa
Shipping	Ice bag