

Kd Validated Anti-Rab27A Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P51159

Item Number: CM0006022



Introduction

KD validated anti-RAB27A rabbit monoclonal antibody for Western blotting (WB), immunofluorescence (IF), immunohistochemistry (IHC), and ELISA. Recognizes human, mouse, and rat RAB27A. Ideal for membrane trafficking and exocytosis studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] RAB27A Rabbit mAb
Remarks/Alias	GS2; RAM; RAB27; HsT18676; RAB27A
Species	Human
GeneID (Human)	5873
GeneID	5873
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human RAB27A (NP_004571.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51159
Protein Weight	25kDa
Shipping	Ice bag