


Kd Validated Anti-Rhot1/Miro1 Rabbit Monoclonal Antibody For Wb And Elisa - Q8Ixi2

Item Number: CM0007289



Introduction

KD validated rabbit monoclonal antibody targeting RHOT1/MIRO1, a mitochondrial Rho GTPase involved in trafficking and fission. Suitable for western blot and ELISA on human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] RHOT1/MIRO1 Rabbit mAb
Remarks/Alias	ARHT1; MIRO1; MIRO-1; [KO/KD Validated] RHOT1/MIRO1
Species	Human
GeneID (human)	55288
GeneID	55288
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-590 of human RHOT1/MIRO1 (NP_060777.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IXI2
Protein Weight	71kDa
Shipping	Ice bag