


Kd Validated Anti-Rhot1/Miro1 Rabbit Polyclonal Antibody - Q8Ixi2

Item Number: CM0024180



Introduction

KD validated rabbit polyclonal antibody targeting human RHOT1/MIRO1, an atypical mitochondrial NTPase. Applications: WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Molecular weight: 71kDa. Useful for mitochondrial trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] RHOT1/MIRO1 Rabbit pAb
Remarks/Alias	ARHT1; MIRO1; MIRO-1; [KD Validated] RHOT1/MIRO1
Species	Human
GeneID (human)	55288
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 310-590 of human RHOT1/MIRO1. (NP_060777.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IXI2
Protein Weight	71kDa
Shipping	Ice bag