


[Kd Validated] Anti-Serpina3 Rabbit Polyclonal Antibody For Wb, Elisa - P01011

Item Number: CM0028229



Introduction

Rabbit polyclonal antibody against human SERPINA3 (Alpha-1-antichymotrypsin), validated for Western blot and ELISA. Suitable for studying acute phase responses and protease inhibition.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SERPINA3 Rabbit pAb
Remarks/Alias	ACT; AACT; GIG24; GIG25; Alpha-1 Antichymotrypsin (SERPINA3)
Species	Human
Gene ID (Human)	12
Gene ID	12
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01011
Protein Weight	48kDa
Shipping	Ice bag