


Kd Validated Anti-Sharpin Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9H0F6

Item Number: CM0006685



Introduction

KD validated rabbit monoclonal antibody against human Sharpin (SIPL1). Suitable for WB, IF/ICC, ELISA. Cross-reacts with human and rat. Useful for studying NF-kB signaling, inflammation, and linear ubiquitination. SWISS: Q9H0F6.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Sharpin Rabbit mAb
Remarks/Alias	SIPL1; [KD Validated] Sharpin
Species	Human
GeneID (Human)	81858
GeneID	81858
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-126 of human SHARPIN (NP_112236.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H0F6
Protein Weight	40kDa
Shipping	Ice bag