


[Kd Validated] Anti-Sod1 Rabbit Polymab® For Wb, Ip, If, Ihc, Elisa - P00441

Item Number: CM0006091



Introduction

High-quality rabbit monoclonal antibody against human SOD1 (Superoxide dismutase [Cu-Zn]), KD validated. Reacts with human, mouse, and rat. Suitable for Western blot, IP, IF/ICC, IF-P, IHC-P, and ELISA applications in oxidative stress and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SOD1 Rabbit PolymAb®
Remarks/Alias	ALS; SOD; ALS1; IPOA; STAHP; hSod1; HEL-S-44; homodimer
Species	Human
GeneID (human)	6647
GeneID	6647
Immunogen	Recombinant protein of human SOD1; Recombinant fusion protein containing a sequence corresponding to amino acids 2-154 of human SOD1 (NP_000445.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00441
Protein Weight	16kDa
Shipping	Ice bag