

[Kd Validated] Anti-Srebf1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Wtn3

Item Number: CM0012129

Introduction

KD validated rabbit monoclonal antibody targeting SREBF1, suitable for WB, IHC-P, ELISA. Recognizes human and mouse SREBF1, a key transcription factor regulating cholesterol biosynthesis and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SREBF1 Rabbit mAb
Remarks/Alias	ADD1; SREBP1; bHLHd1
Species	Mouse
Gene ID (Human)	6720
Gene ID (Mouse)	20787
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse SREBP1 (NP_035610.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9WTN3
Protein Weight	120 kDa
Shipping	Ice bag