

Kd Validated Anti-Sult2B1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O00204

Item Number: CM0006023



Introduction

KD Validated rabbit monoclonal antibody against human sulfotransferase 2B1 (SULT2B1). Suitable for WB, IHC-P, and ELISA. SULT2B1 is key in steroid and cholesterol metabolism, with roles in skin differentiation and disease. Ideal for metabolism and dermatology research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SULT2B1 Rabbit mAb
Remarks/Alias	HSST2; ARCI14
Species	Human
GeneID (human)	6820
GeneID	6820
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-365 of human SULT2B1 (NP_814444.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O00204
Protein Weight	41kDa
Shipping	Ice bag