


Kd Validated Anti-Tfeb Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa, Chip - P19484

Item Number: CM0008749



Introduction

KD Validated rabbit monoclonal antibody against human TFEB (Transcription factor EB). Detects endogenous levels of TFEB. Cross-reacts with mouse and rat. Suitable for Western blot, IHC-P, ELISA, and ChIP applications. Ideal for autophagy and lysosome research.

[Learn More](#)

Parameter	Value
Product Name	TFEB Rabbit mAb [KD Validated]
Remarks/Alias	TCFEB; BHLHE35; ALPHATFEB
Species	Human
Gene ID (Human)	7942
Gene ID	7942
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-173 of human TFEB (NP_009093.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P19484
Protein Weight	44kDa/53kDa
Shipping	Ice bag