

Kd Validated Anti-Tfeb Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P19484

Item Number: CM0025516



Introduction

Rabbit polyclonal antibody against human TFEB, validated for knockdown specificity. Suitable for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Ideal for autophagy and lysosomal biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] TFEB Rabbit pAb
Remarks/Alias	TCFEB; BHLHE35; ALPHATFEB; TFEB
Species	Human
GeneID (Human)	7942
GeneID	7942
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-173 of human TFEB (NP_009093.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19484
Protein Weight	53 kDa
Shipping	Ice bag