


[Kd Validated] Anti-Timm10 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62072

Item Number: CM0005575



Introduction

Rabbit monoclonal antibody against human TIMM10 (TIM10), a mitochondrial import chaperone. Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying mitochondrial inner membrane translocase and related pathways.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] TIMM10 Rabbit mAb
Remarks/Alias	TIM10; TIM10A; TIMM10A; [KD Validated] TIMM10
Species	Human
GeneID (Human)	26519
GeneID	26519
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-90aa of human TIMM10 (NP_000204.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62072
Protein Weight	10kDa
Shipping	Ice bag