

[Kd Validated] Anti-Tradd Monoclonal Antibody For Wb, Ip, Elisa - Q15628

Item Number: CM0011574



Introduction

KD Validated rabbit monoclonal antibody against human TRADD (TNFRSF1A-associated via death domain). Validated for WB, IP, ELISA; reactive with human, mouse, and rat. Ideal for apoptosis and NF-κB signaling studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] TRADD Rabbit mAb
Remarks/Alias	Hs.89862
Species	Human
GeneID (Human)	8717
GeneID	8717
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 199-312 of human TRADD (NP_003780.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15628
Protein Weight	34 kDa
Shipping	Ice bag