


## Kd Validated Anti-Vimentin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08670

Item Number: CM0017571



### Introduction

High-quality KD validated rabbit polyclonal antibody targeting human Vimentin (P08670). Suitable for WB, IHC-P, IF/ICC, and ELISA applications with cross-reactivity to mouse and rat. Ideal for cytoskeleton and cell migration research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] Vimentin Rabbit pAb                                                                                  |
| Remarks/Alias    | CTRCT30; HEL113; Vimentin; VIM; vimentin; in                                                                        |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 7431                                                                                                                |
| Gene ID          | 7431                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human Vimentin (NP_003371.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | P08670                                                                                                              |
| Protein Weight   | 54kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |