

## [Kd Validated] Anti-Xpnpep3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Nqh7

Item Number: CM0011692



### Introduction

[KD Validated] rabbit monoclonal antibody against XPNPEP3 (APP3/ICP55/NPHPL1), reacts with Human, Mouse, Rat. Suitable for western blot, IHC-Paraffin, and ELISA. Target is a mitochondrial aminopeptidase involved in peptide processing and TNFRSF1B signaling; relevant for nephronophthisis studies.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] XPNPEP3 Rabbit mAb                                                                                   |
| Remarks/Alias    | APP3; ICP55; NPHPL1                                                                                                 |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 63929                                                                                                               |
| GeneID           | 63929                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 32-230 of human XPNPEP3 (NP_071381.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q9NQH7                                                                                                              |
| Protein Weight   | 57kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |