


[Kd Validated] Anti-Ythdc1 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q96Mu7

Item Number: CM0006052



Introduction

Rabbit monoclonal antibody against human YTHDC1 (KD validated), suitable for WB, IHC-P, IP, and ELISA. Cross-reacts with mouse and rat. Detects YTHDC1 at ~85 kDa, a key m6A reader and splicing regulator. Ideal for RNA modification and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] YTHDC1 Rabbit mAb
Remarks/Alias	YT521; YT521-B; [KD Validated] YTHDC1
Species	Human
Gene ID (Human)	91746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-71 of human YTHDC1 (NP_001026902.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96MU7
Protein Weight	85kDa (observed)
Shipping	Ice bag