


Kd Validated Anti-Ythdc2 Rabbit Polyclonal Antibody For Wb, Ip, If-P, Elisa - Q9H6S0

Item Number: CM0021397



Introduction

KD validated rabbit polyclonal antibody targeting human YTHDC2 (3'-5' RNA helicase) for WB, IP, IF-P, ELISA. Cross-reacts with mouse and rat. Detects ~160 kDa band; ideal for meiosis and m6A RNA decay research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] YTHDC2 Rabbit pAb
Remarks/Alias	CAHL; hYTHDC2; YTHDC2
Species	Human
GeneID (Human)	64848
GeneID	64848
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1286-1430 of human YTHDC2 (NP_073739.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H6S0
Protein Weight	160kDa
Shipping	Ice bag