

[Kd Validated] Anti-Znf143 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, And Elisa - P52747

Item Number: CM0007651

Introduction

KD-validated rabbit monoclonal antibody targeting human ZNF143 (SBF/STAF). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for gene regulation and mitochondrial function studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] ZNF143 Rabbit mAb
Remarks/Alias	SBF; STAF; pHZ-1; [KD Validated] ZNF143
Species	Human
Gene ID (Human)	7702
Gene ID	7702
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 592-638 of human ZNF143 (NP_003433.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52747
Protein Weight	69kDa
Shipping	Ice bag