


[Kd Validated] Cdk2 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P24941

Item Number: CM0008540



Introduction

High-quality CDK2 rabbit monoclonal antibody for WB, IF/ICC, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Key for cell cycle, DNA repair, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CDK2 Rabbit mAb
Remarks/Alias	CDKN2; p33(CDK2)
Species	Human
GeneID (Human)	1017
GeneID	1017
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-298 of human CDK2 (NP_001789.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24941
Protein Weight	30 kDa/34 kDa
Shipping	Ice bag