

## [Kd Validated] Cdx2 Rabbit Pab - Q99626

Item Number: CM0021891



### Introduction

CDX2 Rabbit Polyclonal Antibody (KD Validated) for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, rat. Ideal for studying intestinal epithelium transcription factor function.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] CDX2 Rabbit pAb                                                                                   |
| Remarks/Alias    | CDX3; CDX-3; CDX2/AS; CDX2                                                                                       |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 1045                                                                                                             |
| Gene ID          | 1045                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 14-165 of human CDX2 (NP_001256.3) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q99626                                                                                                           |
| Protein Weight   | 34kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |